



# Recombinant Escherichia coli O157:H7 L-seryl-tRNA (Sec) selenium transferase

<b>Product Code</b>	CSB-YP481187EOE
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B5YW91
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7 (strain EC4115 / EHEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTTETRSLYS QLPAIDRLLR DSSFLSLRDT YGHTRVVELL RQMLDEAREV IRGSQTLPAW CENWAQEVDVA RLTKEAQSAL RPVINLTGTV LHTNLGRALQ AEEAVEAVAQ AMRSPVTLEY DLDDAGRGRH DRALAQLLCR ITGAEDACIV NNNAAVLLM LAATASGKEV VVSRGELVEI GGAFRIPDVM RQAGCTLHEV GTTNRTHAND YRQAVNENTA LLMKVHTSNY SIQGFCKAID EAELVALGKE LDVPVVDLG SGSLVDLSQY GLPKEPMPQE LIAAGVSLVS FSGDKLLGGP QAGIIVGKKE MIARLQSHPL KRALRADKMT LAALATLRL YLHPEALSEK LPTLRLLTRS AEVIQIQQR LQAPLAAHYG AEFVQVMPC LSQIGSGSLP VDRLPSAALT FPHDGRGSH LESLAARWRE LPVPVIGRIY DGRLWLDLRC LEDEQRFLEM LLK
<b>Source</b>	Yeast
<b>Target Names</b>	selA
<b>Protein Names</b>	Recommended name: L-seryl-tRNA(Sec) selenium transferase EC= 2.9.1.1 Alternative name(s): Selenocysteine synthase Short name= Sec synthase Selenocysteinyl-tRNA(Sec) synthase
<b>Expression Region</b>	1-463
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.