



Recombinant Escherichia coli O157:H7 Phosphopantetheine adenylyltransferase (coaD)

Product Code	CSB-MP481188EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YWD2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MQKRAIYPGT FDPITNGHID IVTRATQMFD HVILAIASP SKKPMFTLEE RVALAQQATA HLGNEVVGFD SLDMANFARN QHATVLRGL RAVADFEYEM QLAHMNRHLM PELESVFLMP SKEWSFISS LVKEVARHQQ DVTHFLPENV HQALMAKLA
Source	Mammalian cell
Target Names	coaD
Protein Names	Recommended name: Phosphopantetheine adenylyltransferase EC= 2.7.7.3 Alternative name(s): Dephospho-CoA pyrophosphorylase Pantetheine-phosphate adenylyltransferase Short name= PPAT
Expression Region	1-159
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.