



Recombinant Mouse Chymase (Cma1), partial

Product Code	CSB-BP005599MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21844
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	IIGGTECIP HSRPYMAYLE IVTSENYLSA CSGFLIRRNF VLTAAHCAGR SITVLLGAHN KTSKEDTWQK LEVEKQFLHP KYDENLVVHD IMLLKLKEKA KLTLGVGTLP LSAFNFP FIP GRMCRAVGWG RTNVNPPASD TLQEVKMLRQ EPQACKHFTS FRHNSQLCVG NPKKMQRVYK GDSGGPLLCA GIAQGIASYV HRNAKPPAVF TRISHYRPWI NKILRE
Source	Baculovirus
Target Names	Cma1
Protein Names	Recommended name: Chymase EC= 3.4.21.39 Alternative name(s): Alpha-chymase Mast cell chymase 1 Mast cell protease 5 Short name= mMCP-5 Mast cell protease I
Expression Region	22-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.