



Recombinant Escherichia coli 3-isopropylmalate dehydratase large subunit (leuC)

Product Code	CSB-MP482307ENM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7L4J3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain 55989 / EAEC)
Purity	>85% (SDS-PAGE)
Sequence	MAKTLYEKLF DAHVVYEAEN ETPLLYIDRH LVHEVTSPQA FDGLRTHGRP VRQPGKTFAT MDHNVSTQTK DINACGEMAR IQMQELIKNC KEFGVELYDL NHPYQGIVHV MGPEQGVTLT GMTIVCGDSH TATHGAFGAL AFGIGTSEVE HVLATQTLKQ GRAKTMKIEV QGKAAPGITA KDIVLAIIGK TGSAGGTGHV VEFCGEAIRD LSMEGRMTLC NMAIEMGAKA GLVAPDETTF NYVKGRLHAP KKGDFDDAVA YWKTQLTDEG ATFDTVVTLQ AEEISPVQVW GTNPGQVISV NDNIPDPASF ADPVERASAE KALAYMGLKP GIPLTEVAID KVFIGSCTNS RIEDLRAAAE IAKGRKVAPG VQALVVPVSG PVKAQAEAEAG LDKIFIEAGF EWRLPGCSMC LAMNNDRLNP GERCASTSNR NFEGRQGRGG RTHLVSPAMA AAAAVTGHFA DIRNIK
Source	Mammalian cell
Target Names	leuC
Protein Names	Recommended name: 3-isopropylmalate dehydratase large subunit EC= 4.2.1.33 Alternative name(s): Alpha-IPM isomerase Short name= IPMI Isopropylmalate isomerase
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.