



# Recombinant Escherichia fergusonii Guanosine-5'-triphosphate,3'-diphosphate pyrophosphatase (gppA)

<b>Product Code</b>	CSB-EP482469EOR-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7LU76
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGSTSSLYAA IDLGSNSFHM LVVREVAGSI QTLTRIKRKV RLAAGLNSEN ALSNEAMERG WQCLRLFAER LQDIPPSQIR VVATATLRLA VNAGDFIAKA QEILGCPVQV ISGEEEARLI YQGVAHTTGG ADQRLVVDIG GASTELVTGT GAQTTSLSL SMGCVTWLER YFADRNLGQE NFDAAEKAAR EVLRPVADEL RYHGWKVCVG ASGTVQALQE IMMAQGMDER ITLEKLQQLK QRAIHCGRLE ELEIDGLTLE RALVFPSGLA ILIAIFTELN IQCMTLAGGA LREGLVYGML HLAVEQDIRS RTLRNIQRRF MIDIDQAQRV AKVAANFFDQ VENEWHLEAI SRDLLISACQ LHEIGLSVDF KQAPQHAAYL VRNLDLPGFT PAQKLLATL LLNQTNPIDL SSLHQQNAVPRV AEQLCRL LRLAIFASR RRDDLVPENT LKANHELLTL TLPQGWLTH PLGKEIIDQE SQWQSYVHWP LEVH
<b>Source</b>	E.coli
<b>Target Names</b>	gppA
<b>Protein Names</b>	Recommended name: Guanosine-5'-triphosphate,3'-diphosphate pyrophosphatase EC= 3.6.1.40 Alternative name(s): Guanosine pentaphosphate phosphohydrolase pppGpp-5'-phosphohydrolase
<b>Expression Region</b>	1-494
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.