



Recombinant Escherichia coli O81 Electron transport complex protein RnfC (rnfC)

Product Code	CSB-YP482737EOP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7MVA7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O81 (strain ED1a)
Purity	>85% (SDS-PAGE)
Sequence	MLKLFSAFRK NKIWDFNGGI HPPEMKTQSN GTPLRQVPLA QRFVIPLKQH IGAEGELCVS VGDKVLRGQP LTRGRGKMLP VHAPTSGLVT AIAPHSTAHP SALAELSVII DADGEDCWIP RDGWADYRSR SREELIERIH QFGVAGLGGA GFPTGVKLQG GGDKIETLII NAAECEPYIT ADDRMLQDCA AQVVEGIRIL AHILQPREIL IGIEDNKPQA ISMLRAVLAD SHDISLRVIP TKYPSGGAKQ LTYILTQKQV PHGGRSSDIG VLMQNVGTAY AVKRAVIDGE PITERVVTLT GEAIARPGNV WARLGPVVRH LLNDAGFCPS ADQMVMGGP LMGFTLPWLD VPVVKITNCL LAPSANELGE PQEEQSCIRC SACADACPAD LLPQQLYWFS KGQQHDKATT HNIADCIECG ACAWVCPSNI PLVQYFRQEK AEIAAIRQEE KRAAEAKARF EARQARLERE KAARLERHKS AAVQPAKDK DAIAAALARV KEKQAQATQP IVIKAGERPD NSAIIAAREA RKAQARAKQA ELQQTNDAAT VADPRKTAVE AAIARAKARK LEQQQANAEP EEQIDPRKAA VEAIIARAKA RKLEQQQANA EPEEQIDPRK AAVEAAIARA KARKLEQQQA NAEPEEQIDP RKAAVEAAIA RAKARKLEQQ QANAEPPEEQI DPRKAAVAAA IARVQAKKAA QQKVVNED
Source	Yeast
Target Names	rsxC
Protein Names	Recommended name: Electron transport complex protein RnfC
Expression Region	1-708
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.