



# Recombinant Escherichia coli O81 DNA ligase (ligA)

<b>Product Code</b>	CSB-EP482754EOP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7MY64
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O81 (strain ED1a)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESIEQQLTE LRTTLRHHEY LYHVMDAPEI PDAEYDRLMR ELRELETKHP ELITPDSPTQ RVGAAPLAAF SQIRHEVPML SLDNVFDEES FLAFNKRVQD RLKSNEKVTW CCELKLDGLA VSILYENGV L VSAATR GDGT TGEDITSNVR TIRAIPLKLH GENIPARLEV RGEVFLPQAG FEKINEDARR TGGKVFANPR NAAAGSLRQL DPRITAKRPL TFFCYG VGV L EGGELPDTHL GR LMQFKAWG LPVSDRVTLC ESAEEVLAFY HEVEKDRPTL GFDIDGVV I K VNSLAQQEQL GFVARAPRWA VAFKFP AQEQ MTFVRDVEFQ VGRTG AITPV ARLEPVHVAG VLVSNATLHN ADEIERLGLR IGDKVVIRRA GDVIPQVVNV VLSERPEDTR EVVFPTHCPV CGSDVERVEG EAVARCTGGL ICGAQRKESL KH FVSRRAMD VDGMGDKIID QLVEKEYVHT PADLFKLTAG KLTGLERMGL KSAQNVVNAL EKAKETTFAR FLYALGIREV GEATAAGLAA YFGTLEALEA ASIEELQKVP DVGIVVASHV HNFFAEESNR N VISELLAEG VHWPEPIVIN AEEIDSPFAG KTVVLTGSL S QMSRDDAKAR LVELGAKVAG SVSKKTDLVI AGEAAGSKLA KAQELGIEVI DETEMLRLLG S
<b>Source</b>	E.coli
<b>Target Names</b>	ligA
<b>Protein Names</b>	Recommended name: DNA ligase EC= 6.5.1.2 Alternative name(s): Polydeoxyribonucleotide synthase [NAD(+)]
<b>Expression Region</b>	1-671
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.