



Recombinant Escherichia coli O17:K52:H18 Penicillin-insensitive murein endopeptidase (mepA)

Product Code	CSB-EP482804EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7N5U0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	A TPWQKITQPV PGSAQSIGSF SNGCIVGADT LPIQSEHYQV MRTDQRRYFG HPDLVMFIQR LSSQVSNLSL GTVLIGDMGM PAGGRFNGGH ASHQTGLDVD IFLQLPKTRW TSAQLLRPQA LDLVSRDGKH VVPTLWKPEI FSLIKLAAQD KDVTRIFVNP AIKQQLCLDA GTDRDWLRKV RPWFQHRAHM HVRLRCPADS LECEDQPLPP PGDGCGAELQ SWFEPKPGT TKPEKTPPP LPPSCQALLD EHVI
Source	E.coli
Target Names	mepA
Protein Names	Recommended name: Penicillin-insensitive murein endopeptidase EC= 3.4.24.- Alternative name(s): D-alanyl-D-alanine-endopeptidase Short name= DD- endopeptidase
Expression Region	20-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.