



# Recombinant Escherichia coli O17:K52:H18 3- (3-hydroxy-phenyl)propionate/3-hydroxycinnamic acid hydroxylase

<b>Product Code</b>	CSB-MP482831EOG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7N8Q4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAIQHPDIQP AVNHSVQVAI AGAGPVGLMM ANYLGQMGID VLVVEKLDKL IDYPRAIGID DEALRTMQSV GLVDNVLPH TTPWHAMRFLT PKGRCFADIQ PMTDEFGWPR RNAFIQPQVD AVMLEGLSRF PNVRCFLFSRE LEAFSQQDDE VTLHLKTEEG QRETVKAQWL VACDGGASFV RRTLNVPFEG KTAPNQWIVV DIANDPLSTP HIYLCCDPVR PYVSAALPHA VRRFEFMVMP GETEEQLREP QNMRKLLSKV LPNPDNVELI RQRVYTHNAR LAQRFRIDRV LLAGDAAHIM PVWQQGGYNS GMRDAFNLA W KLALVIQGKA RDALLDTYQQ ERRDHAKAMI DLSVTAGNVL APPKRWQGT L RDGVSWLLNY LPPVKRYFLE MRFKPM PQYY GGALVREGEA KHSPVGK MFI QPKVTLENGD VTL LDNAIGA NFAVIGWG CN PLWGM SDEQI QQWRALGTRF IQVVPEVQIH TAQDNHDGVL RVGDTQGR LR SWFAQHNASL VVMRPDRFVA ATAIPQTLGK TLNKLASVMT LTRPADADSV EKVA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	mhpA
<b>Protein Names</b>	Recommended name: 3-(3-hydroxy-phenyl)propionate/3-hydroxycinnamic acid hydroxylase Short name= 3-HCI hydroxylase Short name= 3-HPP hydroxylase EC= 1.14.13.127
<b>Expression Region</b>	1-554
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.