



# Recombinant Escherichia coli O17:K52:H18 Ketol-acid reductoisomerase (ilvC)

<b>Product Code</b>	CSB-YP482875EOG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7NF81
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MANYFNTLNL RQQLAQLGKC RFMGRDEFAD GASYLQGKKV VIVGCGAQQGL NQGLNMRDSG LDISYALRKE AIAEKRASWR KATENGFKVG TYEELIPQAD LVVNLTPDKQ HSDVVRTVQP LMKDGAALGY SHGFNIVEVG EQIRKDITVV MVAPKCPGTE VREEYKRGFG VPTLIAVHPE NDPKGEGMAI AKAWAAATGG HRAGVLESSF VAEVKSDLMG EQTILCGMLQ AGSLLCFDKL VEEGTDPAYA EKLIQFGWET ITEALKQGGI TLMMDRLSNP AKLRAYALSE QLKEIMAPLF QKHMDDIISG EFSSGMMADW ANDDKLLTW REETGKTAFE TAPQYEGKIG EQEYFDKGV L MIAMVKAGVE LAFETMVDSG IIEESAYYES LHELPLIANT IARKRLYEMN VVISDTAEYG NYLFSYACVP LLKPFMAELQ PGDLGKAIFE GAVDNAQLRD VNEAIRSHAI EQVGKLRGY MTDMKRIAVA G
<b>Source</b>	Yeast
<b>Target Names</b>	ilvC
<b>Protein Names</b>	Recommended name: Ketol-acid reductoisomerase EC= 1.1.1.86 Alternative name(s): Acetohydroxy-acid isomeroeductase Alpha-keto-beta-hydroxylacil reductoisomerase
<b>Expression Region</b>	1-491
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.