



# Recombinant Escherichia coli O7:K1 3-dehydroquinate synthase (aroB)

<b>Product Code</b>	CSB-EP482928EON
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7NMF2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MERIVVTLGE RSYPTIASG LFNEPASFLP LKSGEQVMLV TNETLAPLYL DKVRGVLEQA GNVVDSVILP DGEQYKSLAV LDTVFTALLQ KPHGRDITTLV ALGGGVVVDL TGFAAASYQR GVRFIQVPTT LLSQVDSSVG GKTAVNHPLG KNMIGAFYQP ASVVVDLDCL KTLPPRELAS GLAEVIKYGI ILDGAFFNWL EENLDALLRL DGPAMAYCIR RCCELKAEVV AADERETGLR ALLNLGHTFG HAIEAEMGYG NWLHGEAVAA GMVMAARTSE RLGQFSSAET QRIITLLTRA GLPVNGPREM SAQAYLPHML RDKKVLAGEM RLILPLAIGK SEVRSGVSHE LVLNAIADCQ SA
<b>Source</b>	E.coli
<b>Target Names</b>	aroB
<b>Protein Names</b>	Recommended name: 3-dehydroquinate synthase EC= 4.2.3.4
<b>Expression Region</b>	1-362
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.