



Recombinant Escherichia coli O7:K1 Glucosamine-6-phosphate deaminase (nagB)

Product Code	CSB-MP482930EON
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NMM9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MRLIPLTTAE QVGKWAARHI VNRINAFKPT ADRPFVLGLP TGGTPMPTYK ALVEMHKAGQ VSFKHVVTFN MDEYVGLPKE HPESYYSFMH RNFFDHVDIP AENINLLNGN APDIDAECRQ YEEKIRSYGK IHLFMGGVGN DGHI AFNEPA SSLASRTRIK TLTHDTRVAN SRFFDNDVNQ VPKYALTVGV GTLLDAEEVM ILVLGSQKAL ALQAAVEGCV NHMWTISCLQ LHPKAIMVCD EPSTMELKVK TLRYFNELEA ENIKGL
Source	Mammalian cell
Target Names	nagB
Protein Names	Recommended name: Glucosamine-6-phosphate deaminase EC= 3.5.99.6 Alternative name(s): GlcN6P deaminase Short name= GNPDA Glucosamine-6-phosphate isomerase
Expression Region	1-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.