



Recombinant Escherichia coli O7:K1 3-demethylubiquinone-9 3-methyltransferase (ubiG)

| | |
|--------------------------|---|
| Product Code | CSB-EP482933EON-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B7NN47 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O7:K1 (strain IAI39 / ExPEC) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MNTENLLKNN NVDHEEIAKF EAVASRWDL EGEFKPLHRI NPLRLGYIAE RAGGLFGKKV LDVGC GGIL AESMAREGAT VTGLDMGFEP LQVAKLHALE SGIQVEYVQE TVEEHAAKHA GQYDVVTCME MLEHVPDPQS VVRACAQLVK PGGDVFFSTL NRNGKSWLMA VVGAEYILRM VPKGTHDVKK FIKPAELLGW VDQTSKERH MTGLHYNPIT NSFCLGPGVD VNYMLHTQNK |
| Source | E.coli |
| Target Names | ubiG |
| Protein Names | Recommended name: 3-demethylubiquinone-9 3-methyltransferase EC= 2.1.1.64 Alternative name(s): 2-polyprenyl-6-hydroxyphenyl methylase EC= 2.1.1.222 3,4-dihydroxy-5-hexaprenylbenzoate methyltransferase Short name= DHMB meth |
| Expression Region | 1-240 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |