



Recombinant Escherichia coli O127:H6 tRNA1 (Val) (adenine (37)-N6)-methyltransferase

| | |
|--------------------------|---|
| Product Code | CSB-EP483003EOB-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B7UH15 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O127:H6 (strain E2348/69 / EPEC) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSQSTSVFRR NGFTFKQFFV AHDRCAMKVG TDGILLGAWA PVAGVKRCLD IGAGSGLLAL MLAQRDDSV MIDAVELESE AAAQAQENIN QSPWAERINV HTADILQWIT QQTVRFDLII SNPPYYQQGV ECATPQREQA RYTTTLDHPS LLTCAAECIT EEGFFCVVLP EQIGNGFTEL ALSMGWHLRL RTDVAENEAR LPHRVLLAFS PQAGECFSDR LVIRGPDQNY SEAYTALTQA FYLFM |
| Source | E.coli |
| Target Names | yfiC |
| Protein Names | Recommended name: tRNA1(Val) (adenine(37)-N6)-methyltransferase EC= 2.1.1.223 Alternative name(s): tRNA m6A37 methyltransferase |
| Expression Region | 1-245 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |