



Recombinant Escherichia coli O127:H6 Glutamyl-tRNA reductase (hemA)

Product Code	CSB-EP483077EOB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7UQ98
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTLLALGINH KTAPVSLRER VSFSPDKLDQ ALDSLLAQPM VQGGVVLSTC NRTELYLSVE ERDDLQEALI RWLCDYHNLN EDDLRLNSLYW HQDNDAVSHL MRVASGLDSL VLGEPQILGQ VKKAFADSQK GHMKASELER MFQKSFSVAK RVRTETDIGA SAVSVAFAAC TLARQIFESL STVTVLLVGA GETIELVARH LREHKVQKMI IANRTRERAQ ILADEVGAEV IALSDIDERL READIIISST ASPLPIIGKG MVERALKSRR NQPMLLDIA VPRDVEPEVG KLANAYLYSV DDLQSIISHN LAQRKAAAVE AETIVAQEAS EFMAWLRAQS ASETIRDYRS QAEHVRDELT AKALAALEQG GDAQAIMQDL AWKLTNRLIH APTKSLQQAA RDGDNERLNI LRDSLGL
Source	E.coli
Target Names	hemA
Protein Names	Recommended name: Glutamyl-tRNA reductase Short name= GluTR EC= 1.2.1.70
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.