



Recombinant Escherichia fergusonii 4-hydroxy-3-methylbut-2-en-1-yl diphosphate synthase (ispG)

Product Code	CSB-MP483268EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LKC3
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MHNQAPIQRR KSTRIYVGNV PIGDGAPIAV QSMTNTRTTD VEATVNQIKA LERVGADIVR VSVPTMDAAE AFKLIKQQVN VPLVADIHFD YRIALKVAEY GVDCLRINPG NIGNEERIRM VVDCARDKNI PIRIGVNAGS LEKDLQEKYG EPTPQALLES AMRHVDHLDR LNFDQFKVSV KASDVFLAVE SYRLAKQID QPLHLGITEA GGARSGAVKS AIGLGLLLSE GIGDTRLVSL AADPVVEEIKV GFDILKSLRI RSRGINFIAC PTCSRQFDV IGTVNALEQR LEDIITPMDV SIIGCVVNGP GEALVSTLGV TGGNKKSGLY EDGVRKDRLD NNDMIDQLEA RIRAKASQLD EARRIDVQQV EK
Source	Mammalian cell
Target Names	ispG
Protein Names	Recommended name: 4-hydroxy-3-methylbut-2-en-1-yl diphosphate synthase EC= 1.17.7.1 Alternative name(s): 1-hydroxy-2-methyl-2-(E)-butenyl 4-diphosphate synthase
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.