



Recombinant Escherichia fergusonii Quinolinate synthase A (nadA)

Product Code	CSB-EP483271EOR-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LKL5
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	≥85% (SDS-PAGE)
Sequence	MSVMFDPDTA IYPFPTKPTP LSIDEKAYYR EKIKRLLKER NAVMVAHYTY DPEIQQLAEE TGGCISDSLE MARFGAKHPA STLLVAGVRF MGETAKILSP EKTILMPTLQ AECSLDLGCP VDQFNAFCDA HPDRTVVVYA NTSAVVKARA DWVVTSSIAV ELIDHLDLGLG EKIIWAPDKH LGRYVQKQTG ADILCWQGTG IVHDEFKTQA LTRLQKQYPD AAILVHPESP QAIVDMADAV GSTSQLIAAA KTLPHQRLIV ATDRGIFYKM QQAVPDKELL EAPTAGEGAT CRSCAHC PWM AMNGLQAI AE ALEQEE SNYE VHVDERLRER ALVPLNRMLD FAATLRG
Source	E.coli
Target Names	nadA
Protein Names	Recommended name: Quinolinate synthase A EC= 2.5.1.72
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.