



# Recombinant Escherichia fergusonii Ribosomal RNA small subunit methyltransferase B (rsmB)

<b>Product Code</b>	CSB-MP483304EOR
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7LRQ5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKKQRNLRSM AAQAVEQVVE QGQSLSNILP PLQQKVSDKD KALLQELCFG VLRTLSQLDW LINKLMARPM TGKQRTVHYL IMVGLYQLLY TRIPPHAALA ETVEGAVAIK RPQLKGLING VLRQFQRQQE ELLAEFNASD ARYLHPSWLL KRLQKAYPEQ WQSIIEANNQ RPPMWLRVNR THHSRDSWLA LLDEAGMKGF PHADYPAVR LETPAPVHVL PGFEEGWVTV QDASAQGCMA WLAPQNGEHI LDLCAAPGGK TTHILEVAPE AQVVAVDIDE QRLSRVYDNL KRLGMKATVK QGDGRYPSQW CGEQQFDRIL LDAPCSATGV IRRHPDIKWL RRDRDIPELA QLQSEILDAI WPHLKSGGIL IYATCSVLPE ENSLQIKAFI QRTADAELCE TGTPEQPGKQ NLPGAEEDG FFYAKLIKK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rsmB
<b>Protein Names</b>	Recommended name: Ribosomal RNA small subunit methyltransferase B EC= 2.1.1.176 Alternative name(s): 16S rRNA m5C967 methyltransferase rRNA (cytosine-C(5)-)-methyltransferase RsmB
<b>Expression Region</b>	1-429
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.