



# Recombinant Escherichia fergusonii RNA 3'-terminal phosphate cyclase (rtcA)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP483312EOR-B   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B7LSC2  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MKRMIALDGA QGEGGGQILR SALSLSMITG QPFTITGIRA GRAKPGLLRQ<br>HLTAVKAATE ICGATVEGAE LGSQRLVFRP GTVRGGDYRF AIGSAGSCTL<br>VLQTVLPALW FADGPSRVEV SGGTDNPSAP PADFIRRVLE PLLAKIGIHQ<br>QTTLLRHGFY PAGGGVVATE VSPVASFNTL QLGERGNIVQ MRGEVLLAGV<br>PRHVAEREIA TLAGSFSLHE QNIHNLPRDQ GPGNTVSLEV ESENIERFF<br>VVGEKRVSAAE VVAAQLVQEV KRYLASPAAV GEYLADQLVL PMALVGAGEF<br>TVAHPSCHLL TNIAVVERFL QVRFCLIEDT GVTRVSIK       |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | rtcA  |
| <b>Protein Names</b>     | Recommended name: RNA 3'-terminal phosphate cyclase Short name= RNA cyclase Short name= RNA-3'-phosphate cyclase EC= 6.5.1.4  |
| <b>Expression Region</b> | 1-338   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |