



Recombinant Escherichia fergusonii UPF0229 protein yeaH (yeaH)

Product Code	CSB-MP483316EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LSN5
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MTWFIIDRRLLN GKNKSMVNRQ RFLRRYKAQI KQSISEAINK RSVTDVDSGE SVSIPTEDIS EPMFHQGRGG LRHRVHPGND HFVQSDRIER PQGGGGGSGS GQQQASQDGE GQDEFVVFQIS KDEYLDLLFE DLALPNLKRN QQRQLTEYKT HRAGYTANGV PANISVVRSL QNSLARRTAM TAGKRRELHA LEETLQIISN SEPAQLLEEE RLRKEIAELR EKISRVPFID TFDLRYKNYE KRDPDPSSQAV MFCLMDVSGS MDQSTKDMAK RFYILLYLFL SRTYKNVEVV YIRHHTQAKE VDEHEFFYSQ ETGGTIVSSA LKLMDEVVQA RYDPAQWNIY AAQASDGDNW ADDSPLCHEI LAKKLLPVVR YYSYIEITRR AHQTLWREYE HLQSTFENFA MQHIRDQDDI YPVFRELPHK QNATVKD
Source	Mammalian cell
Target Names	yeaH
Protein Names	Recommended name: UPF0229 protein yeaH
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.