



Recombinant Escherichia fergusonii Phosphate acyltransferase (plsX)

Product Code	CSB-MP483318EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LT48
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MTRLTLALDV MGGDFGPSVT VPAALQALNS NSQLTLLLVG NPDTITPLLA KADFEQRSRL QIIPAQSVIA SNARPSQAIR ASRGSSMRVA LELVKEGRAQ ACVSAGNTGA LMGLAKLLLK PLEGIERPAL VTVLPHQQKG KTVVLDLGAN VDCDSTMLVQ FAIMGSVLAE EVVGIPNPRV ALLNIGEEEV KGLDSIRDAS AVLKTISSVN YIGYLEANEL LTGKTDVLVC DGFTGNVTLK TMEGVVRFML SLLKSQGEGK KRSWWLLLLK RWLQKSLTRR FSHLNPQYN GACLLGLRGT VIKSHGAANQ RAFAVAIEQA VQAVQRQVPQ RIAARLKSIVY PAGFELLEDE RGKKPNHR
Source	Mammalian cell
Target Names	plsX
Protein Names	Recommended name: Phosphate acyltransferase EC= 2.3.1.n2 Alternative name(s): Acyl-ACP phosphotransacylase Acyl-[acyl-carrier-protein]--phosphate acyltransferase Phosphate-acyl-ACP acyltransferase
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.