



# Recombinant Escherichia fergusonii ADP-L-glycero-D-manno-heptose-6-epimerase (hldD)

<b>Product Code</b>	CSB-BP483334EOR
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7LVH8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MIIVTGGAGF IGSNIVKALN DKGITDILVV DNLKDGTKFV NLVDLNIADY            MDKEDFLIQI MAGEEFGDVE AIFHEGACSS TTEWDGKYMM DNNYQYSKEL            LHYCLEREIP FLYASSAATY GGRTSDFIES REYEKPLNVY GYSKFLFDEY            VRQILPEANS QIVGFYFNV YGPREGHKGS MASVAFHLNT QLNNGESPKL            FEGSENFKRD FVYVGDVADV NLWFLENGVS GIFNLGTGRA ESFQAVADAT            LAYHKKGQIE YIPFPDKLKG RYQAFTQADL TNLRAAGYDK PFKTVAEGVT            EYMAWLNLRDA</p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	hldD
<b>Protein Names</b>	Recommended name: ADP-L-glycero-D-manno-heptose-6-epimerase EC=5.1.3.20 Alternative name(s): ADP-L-glycero-beta-D-manno-heptose-6-epimerase Short name= ADP-glyceromanno-heptose 6-epimerase Short name= ADP-hep 6-epimerase
<b>Expression Region</b>	1-310
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.