



# Recombinant Escherichia coli O8 Carboxylate-amine ligase YbdK (ybdK)

<b>Product Code</b>	CSB-EP483405E00
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7M4R1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O8 (strain IAI1)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPLPDFHVSE PFTLGIELEM QVVNPPGYDL SQDSSMLIDA VKNEITAGEV KHDITESMLE LATDVCRDIN QAAGQFSAMQ KVVQLQAAADH HLEICGGGTH PFQKWQRQEV CDNERYQRTL ENFGYLIQQA TVFGQHVHVG CASGDDAIYL LHGLSRFVPH FIALSAASPY MQGTDTRFAS SRPNIFSAPP DNGPMPWVSN WQQFEALFRC LSYTTMIDSI KDLHWDIRPS PHFGTVEVRV MDTPLTLSHA VNMAGLIQAT AHWLLTERPF KHQEKDYLLY KFNRFQACRY GLEGVITDPH TGDRRPLTED TLRLLEKIAP SAHKMGASSA IEALHRQVVS GLNEAQLMRD FVADGGSLIG LVKKHCEIWA GD
<b>Source</b>	E.coli
<b>Target Names</b>	ybdK
<b>Protein Names</b>	Recommended name: Carboxylate-amine ligase YbdK EC= 6.3.-.-
<b>Expression Region</b>	1-372
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.