



# Recombinant Escherichia coli O45:K1 Acetyl-coenzyme A carboxylase carboxyl transferase subunit alpha (accA)

<b>Product Code</b>	CSB-MP483464EOK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7MBG6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O45:K1 (strain S88 / ExPEC)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSLNFLDFEQ PIAELEAKID SLTAVSRQDE KLDINIDEEV HRLREKSVEL TRKIFADLGA WQIAQLARHP QRPYTLDYVR LAFDEFDELA GDRAYADDKA IVGGIARLDG RPVMIIGHQK GREKKEKIRR NFGMPAPEGY RKALRLMQMA ERFKMPIITF IDTPGAYPGV GAEERGQSEA IARNLREMSR LGVPVVCTVI GEGGSGGALA IGVGDKVNML QYSTYSVISP EGCASILWKS ADKAPLAAEA MGIAPRLKE LKLIDSIPE PLGGAHRNPE AMAASLKAQL LTDLADLDVL STEDLKNRRY QRLMSYGYA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	accA
<b>Protein Names</b>	Recommended name: Acetyl-coenzyme A carboxylase carboxyl transferase subunit alpha Short name= ACCase subunit alpha Short name= Acetyl-CoA carboxylase carboxyltransferase subunit alpha EC= 6.4.1.2
<b>Expression Region</b>	1-319
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.