



Recombinant Escherichia coli O45:K1 Protein mtfA (mtfA)

Product Code	CSB-EP483472EOK-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7MCM7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O45:K1 (strain S88 / ExPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MIKWPWKVQE SAHQ TALPWQ EALSIPLLTC LTEQE QSKLV TLAERFLQQK RLVPLQGFEL NSLRSCRIAL LFCLPVLELG LEWLDSFHEV LIYPAPFVVD DEWEDDIGLV HNQRIVQSGQ SWQQGPIVLN WLDIQDSFDA SGFNLIIEHV AHKLDTRNGD RASGVPIPL REVAGWEHDL HAAMNNIQEE IELVGENAAS IDAYAASDPA ECFAVLSEYF FSAPELFAPR FPSLWQRFCQ FYQQDPLQRL HHANDTDSFS ATNVH
Source	E.coli
Target Names	mtfA
Protein Names	Recommended name: Protein mtfA Alternative name(s): Mlc titration factor A
Expression Region	1-265
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.