



# Recombinant Escherichia coli O45:K1 Biotin synthase (bioB)

<b>Product Code</b>	CSB-EP483503EOK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7MGN3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O45:K1 (strain S88 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAHRPRWTLS QVTELFKPL LDLLFEAQQV HRQHFDPRQV QVSTLLSIKT GACPEDCKYC PQSSRYKTGL EAERLMEVEQ VLESARKAKA AGSTRFCMGA AWKNPHERDM PYLEQMVQGV KAMGLEACMT LGTLESQAQ RLANAGLDYY NHNLDTSPF YGNIITRTY QERLDTLEKV REAGIKVCSG GIVGLGETVK DRAGLLLQLA NLPTPPESVP INMLVKVKG T PLADNDDVDA FDFIRTIAVA RIMMPTS YVR LSAGREQMNE QTQAMCFMAG ANSIFYGCKL LTPNPEEDK DLQLFRKLGL NPQQTAVLAG DNEQQQRLEQ ALMTPDTDEY YNAAAL
<b>Source</b>	E.coli
<b>Target Names</b>	bioB
<b>Protein Names</b>	Recommended name: Biotin synthase EC= 2.8.1.6
<b>Expression Region</b>	1-346
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.