



Recombinant Escherichia coli O17:K52:H18 Probable 4-deoxy-4-formamido-L-arabinose- phosphoundecaprenol deformylase ArnD (arnD)

Product Code	CSB-BP483655EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7N5M1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTKVGLRIDV DFRGTREGV PRLLLEILSKH NIQASIFFSV GPDNMGRHLW RLVKPRFLWK MLRSNAASLY GWDILLAGTA WPGKEIGHAN ADIIREAAKH HEVGLHAWDH HAWQAHSGNW DRQTMVDDIA RGLRTLEEII GQPVTCSAAA GWRADQQVIE AKEAFHLRYN SDCRGAMPFR PLLESGKPGT AQIPVTLPTW DEVIGRDVKA EDFNGWLLNR ILRDKGMPVY TIHAEVEGCA YQHNFDLLK RAAQEGVTFC PLSELLSGTL PLGQVVRGNI AGREGWLGCC QIAGSR
Source	Baculovirus
Target Names	arnD
Protein Names	Recommended name: Probable 4-deoxy-4-formamido-L-arabinose- phosphoundecaprenol deformylase ArnD EC= 3.5.1.n3
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.