



# Recombinant Escherichia coli O17:K52:H18 Protein CsiD (csiD)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP483664EOG   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B7N6P3  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MNALTAVQNN AVDSGQDYSYSG FTLIPSAQSP RLELETFTEQ TTKQFLEQVA<br>EWPVQALEYK SFLRFRVGGKI LDDLCANQLQ PLLLKTLLNR AEGALLINAV<br>GVDDVKQADE MVKLATAVAH LIGRSNFDAM SGQYYARFVV KNVDNSDSYL<br>RQPHRVMELH NDGTYYVEEIT DYVLMMKIDE QNMQGGNSLL LHLDDWEHL<br>HYFRHPLARR PMRFAAPPSK NVSKDVFHPV FDVDQQGRP V MRYIDQFVQP<br>KDFEEGVWLS ELSDAIETSK GILSVPVPG KFLINLFW LHGRDRFTPH<br>PDLRRELMRQ RGYFAYATHH YQTHQ                    |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | csiD  |
| <b>Protein Names</b>     | Recommended name: Protein CsiD  |
| <b>Expression Region</b> | 1-325   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |