



# Recombinant Escherichia coli O17:K52:H18 Crotonobetainyl-CoA dehydrogenase (caiA)

<b>Product Code</b>	CSB-YP483673EOG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7N7R4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDFNLNDEQE LfvAgirelM ASENWEAYFA ECDRDSVYPE RFVKALADMG IDSLIPEEH GGLDAGFVTL AAVWMElGRL GAPTIVLYQL PGGFNTFLRE GTQEQIDKIM AFRGTGKQMW NSAITEPGAG SDVGS�KTTY TRRNGKIYLN GSKCFITSSA YTPYIVVMAR DGASPDkPVY TEWFVDMSKP GIKVTKLEKL GLRMDSCCEI TFDDVELDEK DMFGREGNGF NRVKEEFdHE RFLVALTNYG TAMCAFEDAA RYANQRVQFG EAIGRFQliQ EKFAHMAIKL NSMKNMLYEA AWKADNGTIT SGDAAMCKYF CANAAFEVVD SAMQVLGGVG IAGNHRISRF WRDLRVDRVS GGSDemQILT LGRAVLKQYR
<b>Source</b>	Yeast
<b>Target Names</b>	caiA
<b>Protein Names</b>	Recommended name: Crotonobetainyl-CoA dehydrogenase EC= 1.3.99.- Alternative name(s): Crotonobetainyl-CoA reductase
<b>Expression Region</b>	1-380
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.