



Recombinant Escherichia coli O17:K52:H18 33 kDa chaperonin (hslO)

Product Code	CSB-EP483716EOG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NE05
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MPQHDQLHRY LFENFAVRGE LVTVSETLQQ ILENHDYPQP VKNVLAELLV ATSLLTATLK FDGDITVQLQ GDGPMNLAVI NGNNNQMRG VARVQGEIPE NADLKTLVGN GYVVITITPS EGERYQGVVG LEGDTLAACL EDYFMRSEQL PTRLFIRTGD VDGKPAAGGM LLQVMPAQNA QQDDFDHLAT LTETIKTEEL LTLPANEVLW RLYHEEEVTV YDPQDVEFKC TCSRERCADA LKTLPEEVD SILTEDGEID MHCDYCGNH Y LFNAMDIAEI RNNASPADPQ VH
Source	E.coli
Target Names	hslO
Protein Names	Recommended name: 33 kDa chaperonin Alternative name(s): Heat shock protein 33 homolog Short name= HSP33
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.