



Recombinant Escherichia coli O17:K52:H18 3-keto-L-gulonate-6-phosphate decarboxylase UlaD (ulaD)

Product Code	CSB-EP483737EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NGD0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSLPMLQVAL DNQTMDSAYE TTRQIAEEVD IIEVGTILCV GEGVRAVRDL KALYPHKIVL ADAKIADAGK ILSRMCFEAN ADWVTVICCA DINTAKGALD VAKEFNGDVQ IELTGYWTWE QAQQWRDAGI QQVVYHRSRD AQAAGVAWGE ADITAIKRLS DMGFKVTVTG GLALEDLPLF KGIPIHVFA GRSIRDAASP VEAARQFKRS IAELWG
Source	E.coli
Target Names	ulaD
Protein Names	Recommended name: 3-keto-L-gulonate-6-phosphate decarboxylase UlaD EC=4.1.1.85 Alternative name(s): 3-dehydro-L-gulonate-6-phosphate decarboxylase KGPDC L-ascorbate utilization protein D
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.