



Recombinant Human Lactoperoxidase (LPO)

Product Code	CSB-EP013066HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22079
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	TRTAIRNGQVWEESLKRLRQKASLTNVTDPSLDLTSLSLEVCGGAPAPVVRCDCSPYRTITGDCNNRRKPALGAANRALARWLP AEYEDGLSLPFGWTPGKTRNGFPLPLAREVSNKIVGYLNEEGVLDQNRSLLFMQWGQIVDHDLDLDFAPDTELGSSEYKAQCDEYCIQGDNCFPI MFPPNDPKAGTQGKCMPPFFRAGFVCPTPPYKSLAREQINALTSFLDASFVYSSEPSLASRLRNLS SPLGLMAVNQEVS DHGLPYLPYDSKKPSPCEFIN TTARVPCFLAGDSRASEHILLATSH TLFLREHNRLARELKR LN PQWDGEKLYQEARKILGAFVQIITFRDYLPILLGDHMQKWIPPYQGYSESV DPRISNVFTFAFRFGHLEVPSSMFRLDENYQPWGPEPELPLHTLFFNTWRMV KDGGIDPLVRGLLAKKSKLMKQNKMMTGELRNKLFQPTHRIHGFDLAAINTQR CRDHGQPGYNSWRAFCDLSQPQTLEELNTVLKSKMLAKKLLGLYGTPDNIDI WIGAI AEPLVERGRVGPLLACLLGKQFQQIRDGDRFWWENPGVFTNEQKDSL QKMSFSRLVCDNTRITK VPRDPFWANSYPYDFVDCSAIDKLDLSPWASVKN
Source	E.coli
Target Names	LPO
Protein Names	Recommended name: Lactoperoxidase Short name= LPO EC= 1.11.1.7 Alternative name(s): Salivary peroxidase Short name= SPO
Expression Region	81-712
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.