



Recombinant Escherichia coli O7:K1 Rhamnulokinase (rhaB)

Product Code	CSB-BP483830EON
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NUA2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MTFRNCVAVD LGASSGRVML ARYERECRSL TLREIHRFKN GLHSQNGYVT WDVDSLESAI RLGLNKVCEE GIRIDSIGID TWGVDFVLLD LQGKRVGLPV AYRDSRTNGL MTQAQQQLGK RDIYQRSIGI FLPFNTLYQL RALTEQQPEL IPHIAHALLM PDYFSYRLTG KMNWEYTNAT TTQLVNINSD DWDESLLAWS GANKAWFGHP THPGNVIGHW ICPQGNEIPV VAVASHDTAS AVIASPLNGS RAAYLSSGTW SLMGFESQTP FTNDTALAAN ITNEGGAEGR YRVLKNIMGL WLLQRVLQER QINDLPALIA ATQALPACRF TINPNDDRFI NPDEMCSEIQ AACRETAQPI PESDAELARC IFDSLALLYA DVLHELAQLR GEDFSQLHIV GGGCQNALLN QLCADACGIR VIAGPVEAST LGNIGIQLMT LDELNNVDDF RQVVSTTANL TTFTPNPDSE IAHYVAQIHS TRQTKELCA
Source	Baculovirus
Target Names	rhaB
Protein Names	Recommended name: Rhamnulokinase EC= 2.7.1.5 Alternative name(s): Rhamnulose kinase
Expression Region	1-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.