



Recombinant Escherichia coli O127:H6 Deoxyguanosinetriphosphate triphosphohydrolase (dgt)

Product Code	CSB-MP483870EOB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7UIK7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Purity	>85% (SDS-PAGE)
Sequence	<p>MAQIDFRKKI NWHRRYRSPQ GVKTEHEILR IFESDRGR II NSPAIRRLQQ KTQVFPLERN AAVRTRLTHS MEVQQVGRYI AKEILSRLKE LKLEAYGLD ELTGPFESIV EMSCLMHDIG NPPFGHFGEA AINDWFRQRL YPEDAESQPL TDDRCSVAAL RLRDGE EPLN ALRRKIRQDL CHFEGNAQGI RLVHTLMRMN LTWAQVGGIL KYTRPAWWRG ETPETHHYLM KKPGYYLSEE AYIARLRKEL NLALYSRFPL TWIMEAADDI SYCVADLEDA VEKRIFTVEQ LYHHLHEAWG QHEKGSLSFSL VVENAWEKSR SNSLSRSTED QFFMYLRVNT LNKLVPYAAQ RFIDNLPAIF AGTFNHALLE DASECSDLLK LYKNVAVKHV FSHPDVEQLE LQGYRVISGL LEIYRPLLNL SLSDFTELVE KERVKRFP IE TRLFHKLSTR HRLAYVEAVS KLPSDSPEFP LWEYYYRCRL LQDYISGMTD LYAWDEYRRL MAVEQ</p>
Source	Mammalian cell
Target Names	dgt
Protein Names	Recommended name: Deoxyguanosinetriphosphate triphosphohydrolase Short name= dGTP triphosphohydrolase Short name= dGTPase EC= 3.1.5.1
Expression Region	1-505
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.