



Recombinant Escherichia coli O127:H6 Bifunctional protein GImU (glmU)

Product Code	CSB-BP483901EOB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7UMJ5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Purity	>85% (SDS-PAGE)
Sequence	MLNNAMSVVI LAAGKGTRMY SDLPKVLHTL AGKTMVQHVI DAANELGAAH VHLVYGHGGD LLKQALKNDN LNWVLQAEQL GTGHAMQQAA PFFADDEDIL MLYGDVPLIS VETLQRLRDA KPQGGIGLLT VKLDDPTGYG RITRENGKVT GIVEHKDATD EQRQIQEINT GILIANGADM KRWLAKLTNN NAQGEYYITD IALAYQEGR EIVAVHPQRL SEVEGVNRL QLSRLERVYQ SEQAEKLLLA GVMLRDPARF DLRGTLTHGR DVEIDTNVII EGNVTLGHRV KIGAGCVIKN SVIGDDCEIS PYTVVEDANL AAACITGPFA RLRPGAELLE GAHVGNFVEM KKARLGKGSK AGHLTYLGDA EIGDNNVIGA GTITCNYDGA NKFKTIIGDD VFIGSDTQLV APVTVGKGAT IAAGTTVTRN IGENALAIRS VPQAQKEGWR RPVKKK
Source	Baculovirus
Target Names	glmU
Protein Names	Recommended name: Bifunctional protein GImU Including the following 2 domains: UDP-N-acetylglucosamine pyrophosphorylase EC= 2.7.7.23 Alternative name(s): N-acetylglucosamine-1-phosphate uridylyltransferase Glucosamine-1-phosphate
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.