



# Recombinant Escherichia coli O127:H6 Glucans biosynthesis protein D (mdoD)

<b>Product Code</b>	CSB-BP483928EOB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7URH4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ADSDIADG QTQRFDSLIL QSMAHDLAQT AWRGAPRPLP DTLATMTPQA YNSIQYDAEK SLWHNVENRQ LDAQFFHMGM GFRRRVRFMS VDPATHLARE IHFRPELFKY NDAGVDTKQL EGQSDLGFAG FRVFKAPELA RRDVVSFLGA SYFRAVDDTY QYGLSARGLA IDTYTDSKEE FPDFTAFWFD TVKPGATTFT VYALLDSASI TGAYKFTIHC EKNQVIMDVE NHLYARKDIK QLGIAPMTSM FSCGTNERRM CDTIHPQIHD SDRLSMWRGN GEWICRPLNN PQKLQFNAYT DNNPKGFGLL QLDRDFSHYQ DIMGWYNKRP SLWVEPRNKW GKGTIGLMEI PTTGETLDNI VCFWQPEKAV KAGDEFQAFQY RLYWSAQPPV HCPLARVMAT RTGMGGFPEG WAPGEHYPEK WARRFAVDFV GGD LKAAAPK GIEPVITLSS GEAKQIEILY IEPIDGYRIQ FDWYPTSDST DPVDMRMYLR CQGDAISETW LYQYFPPAPD KRQYVDDRVM S
<b>Source</b>	Baculovirus
<b>Target Names</b>	mdoD
<b>Protein Names</b>	Recommended name: Glucans biosynthesis protein D
<b>Expression Region</b>	33-551
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.