



# Recombinant Escherichia fergusonii Anaerobic glycerol-3-phosphate dehydrogenase subunit B (glpB)

<b>Product Code</b>	CSB-EP484168EOR-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7LM86
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRFDTVIMGG GLAGLLCGLQ LQKHGLRCAI VTRGQSALHF SSGSLDLLSH LPDGQAVTDI HSGLESLRQQ APAHPYSLG PQRVLDLACQ AQTLIAESGA QLQGSVEQAH QRITPLGTLR STWLSSPEVP VWPLPAKKIC VVGISGLMDF QAHLAAASLR ELDLSVETAE IELPELDVLR NNATEFRAVN IARFLDNEEN WSLLLDALIP VANTCEMILM PACFGLADDK LWYWLNEKLP CSLMLLPTLP PSVLGIRLQN QLQRQFVRQG GVWMPGDEVK KVTCKNGMVN EIWTRNHADI PLRPRFAVLA SGSFFSGGLV AEREGIREPI VGLDVLQTAT RGEWYKGDFF APQPWQQFGV TTDEALRPSQ AGKTIENLFA IGSVLGGFDP IAQGCGGGVC AVSALHAAQQ IAQRAGGQQ
<b>Source</b>	E.coli
<b>Target Names</b>	glpB
<b>Protein Names</b>	Recommended name: Anaerobic glycerol-3-phosphate dehydrogenase subunit B Short name= Anaerobic G-3-P dehydrogenase subunit B Short name= Anaerobic G3Pdase B EC= 1.1.5.3
<b>Expression Region</b>	1-419
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.