



Recombinant Escherichia fergusonii tRNA (mo5U34)-methyltransferase

Product Code	CSB-YP484179EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LPH3
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MIDFGNFYSL IAKNHLSHWL ETLPAQIASW QREQQHGLFK QWSNAVEFLP EIKPYRLDLL HSVTAESEEP LSAGQIKRIE TLMRNLMPWR KGPFSLYGVN IDTEWRSDWK WDRVLPPLSD LTGRTILDVG CGSGYHMWRM IGAGAHLAGV IDPTQLFLCQ FEAVRKLLGN DQRAHLLPLG IEQLPALKAF DTVFSGMGLV HRRSPLEHLW QLKDQLVNEG ELVLETLID GDENTVLVPG DRYAQMRNVY FIPSALALKN WLKCKGFVDI RIADSVSTTT EEQRRTWMV TESLSDFLDP HDPGKTVEGY PPKRAVLIA RKP
Source	Yeast
Target Names	cmoB
Protein Names	Recommended name: tRNA (mo5U34)-methyltransferase EC= 2.1.1.-
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.