



# Recombinant Escherichia fergusonii Glucose-1-phosphate adenylyltransferase (glgC)

<b>Product Code</b>	CSB-BP484201EOR
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7LSE1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVGLEKNDPL MLARQLPIKS VALILAGGRG TRKDLTNKR AKPAVHFGGK FRIIDFALSN CINSIRRIG VITQYQSHTL VQHIQRGWSF FNEEMNEFVD LLPAQQRMKG ENWYRGTADA VTQNLDIIRR YKAEYVVILA GDHIYKQDYS RMLIDHVEKG ARCTVACMPV PIQEASAFGV MAVDQDEKII EFVEKPANPP SMPNDPTRSL ASMGIVFDA DYLYELLEED DNDENSSHDF GKDIIPKITQ AGMAYAHPPF LSCVQSDPDS EPYWRDVGTL EAYWKANLDL ASVPELDMY DQHWPIRTYN ESLPPAKFVQ DRSGSHGRTL NSLVSGGCVI SGSVVVQSVL FSRVRVNSFC NIDSAVLLPE VVVGSRCLLR RCVIDRACVI PEGMVIKENA EEDARRFYRS EEGIVLVTRE MLRKLGHKQE R
<b>Source</b>	Baculovirus
<b>Target Names</b>	glgC
<b>Protein Names</b>	Recommended name: Glucose-1-phosphate adenylyltransferase EC= 2.7.7.27 Alternative name(s): ADP-glucose pyrophosphorylase Short name= ADPGlc PPase ADP-glucose synthase
<b>Expression Region</b>	1-431
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.