



Recombinant Escherichia coli O8 Phosphate acyltransferase (plsX)

Product Code	CSB-EP484238E00-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LX24
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O8 (strain IA11)
Purity	>85% (SDS-PAGE)
Sequence	MTRLTLALDV MGGDFGPSVT VPAALQALNS NSQLTLLLVG NPDAITPLLA KADFEQRSRL QIIPAQSVIA SDARPSQAIR ASRGSSMRVA LELVKEGRAQ ACVSAGNTGA LMGLAKLLLK PLEGIERPAL VTVLPHQQKG KTVVLDLGAN VDCDSTMLVQ FAIMGSVLAE EVVEIPNPRV ALLNIGEEEV KGLDSIRDAS AVLKTIPSIN YIGYLEANEL LTGKTDVLC DGFTGNVTLK TMEGVVRMFL SLLKSQGEGK KRSWWLLLLK RWLQKSLTRR FSHLNPQYN GACLLGLRGT VIKSHGAANQ RAFAVAIEQA VQAVQRQVPQ RIAARLESVY PAGFELLDGG KSGTLR
Source	E.coli
Target Names	plsX
Protein Names	Recommended name: Phosphate acyltransferase EC= 2.3.1.n2 Alternative name(s): Acyl-ACP phosphotransacylase Acyl-[acyl-carrier-protein]--phosphate acyltransferase Phosphate-acyl-ACP acyltransferase
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.