



# Recombinant Escherichia coli O8 Succinylglutamate desuccinylase (astE)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP484274E00   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B7M1F6  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O8 (strain IA11)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MDNFLALTLT GK KPVITERE INGVRWRWLG DGVLELTPLT PPQGALVISA<br>GIHGNETAPV EMLDALLGAI SHGEIPLRWR LLVILGNPPA LKQGKRYCHS<br>DMNRMFGGRW QLFAESGETC RARELEQCLE DFYDQGKESV<br>RWHLDLHTAI RGS LHPQFGV LPQRDIPWDE KFLTWLGAAG LEALVFHQEP<br>GGTFT HFSAR HFGALACTLE LGKALPFGQN DLRQFAVTAS AIAALLSGES<br>VGIVRTPPLR YRVVSQITRH SPSFEMHMAS DTLNFMPFEK GTLLAQDGEE<br>RFTVTHDVEY VLFPNPLVAL GLRAGLMLEK IS                     |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | astE  |
| <b>Protein Names</b>     | Recommended name: Succinylglutamate desuccinylase EC= 3.5.1.96  |
| <b>Expression Region</b> | 1-322   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |