



# Recombinant Escherichia coli O45:K1 Bis (5'-nucleosyl)-tetrphosphatase, symmetrical

<b>Product Code</b>	CSB-MP484350EOK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7MAH4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O45:K1 (strain S88 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MATYLIGDVH GCYDELIALL HKVEFTPGKD TLWLTGDLVA RGPGLDVLRL YVKSLGDSVR LVLGNHDLHL LAVFAGISRN KPKDRLTPLL EAPDADELLN WLRRQPLLQI DKEKKLVMAH AGITPQWDLQ TAKECARDVE AVLSSDSYYPF FLDAMYGDMP NNWSPELRGL GRLRFITNAF TRMRFCFPNG QLDMYSKESP EEAPAPLKPW FAIPGPVAEE YNIAFGHWAS LEGKGTPEGI YALDTGCCWG GTLTCLRWED KQYFVQPSNR HKDLSEGEAV AS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	apaH
<b>Protein Names</b>	Recommended name: Bis(5'-nucleosyl)-tetrphosphatase, symmetrical EC= 3.6.1.41 Alternative name(s): Ap4A hydrolase Diadenosine 5',5'''-P1,P4-tetrphosphate pyrophosphohydrolase Diadenosine tetrphosphatase
<b>Expression Region</b>	1-282
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.