



# Recombinant Escherichia coli O45:K1 Phosphoglucosamine mutase (glmM)

<b>Product Code</b>	CSB-EP484358EOK-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7MB95
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O45:K1 (strain S88 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSNRKYFGTD GIRGRVGDAP ITPDFVLKLG WAAGKVLARH GSRKIIIGKD TRISGYMLES ALEAGLAAAG LSALFTGPMP TPAVAYLTRT FRAEAGIVIS ASHNPFYDNG IKFFSIDGTK LPDAVEEAIE AEMEKEISCV DSAELGKASR IVDAAGRYIE FCKATFPNEL SLSELKIVVD CANGATYHIA PNVLRELGAN VIAIGCEPNG VNINAEVGAT DVRALQARVL AEKADLGI AF DGDGDRVIMV DHEGNKVDGD QIMYIIAREG LRQGQLRGGV VGTLMNSMGL ELALKQLGIP FARAKVGDY VLEKMKEKW RIGAENSGHV ILLDKTTTGD GIVAGLQVLA AMARNHMSLH DLCSGMKMFP QILVNVRYTA GSGDPLEHES VKAVTAEVEA ALGSRGRVLL RKSGTEPLIR VMVEGEDEAQ VTEFAHRIAD AVKAV
<b>Source</b>	E.coli
<b>Target Names</b>	glmM
<b>Protein Names</b>	Recommended name: Phosphoglucosamine mutase EC= 5.4.2.10
<b>Expression Region</b>	1-445
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.