



Recombinant Escherichia coli O45:K1 Sulfate adenylyltransferase subunit 2 (cysD)

Product Code	CSB-MP484431EOK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7MKM6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O45:K1 (strain S88 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MDQKRLTHLR QLEAESIHII REVAAEFSNP VMLYSIGKDS SVMLHLARKA FYPGTLPPFPL LHVDTGWKFR EMYEFRDRTA KAYGCELLVH KNPEGVAMGI NPFVHGSAXH TDIMKTEGLK QALNKYGFDA AFGGARRDEE KSKRAKERIYS FRDRFHRWDP KNQRPELWHN YNGQINKGES IRVFPLSNWT EQDIWQYIWL ENIDIVPLYL AAERPVLERD GMLMMIDDNR INLQSGEVIK KRMVRFRTL G CWPLTGAVES NAQTLPEIIE EMLVSTTSER QGRVIDRDQA GSMELKKRQG YF
Source	Mammalian cell
Target Names	cysD
Protein Names	Recommended name: Sulfate adenylyltransferase subunit 2 EC= 2.7.7.4 Alternative name(s): ATP-sulfurylase small subunit Sulfate adenylate transferase Short name= SAT
Expression Region	1-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.