



# Recombinant Escherichia coli O17:K52:H18 Arginine N-succinyltransferase (astA)

<b>Product Code</b>	CSB-BP484552EOG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7N583
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMVIRPVERS DVSALMQLAS KTGGLTSLP ANEATLSARI ERAIKTWQGE LPKSEQGYVF VLEDSETGTV AGICAIEVAV GLNDPWYNYR VGTLVHASKE LNVYNALPTL FLSNDHTGSS ELCTLFLDPE WRKEGNGYLL SKSRFMFMAA FRDKFNDKVV AEMRGVIDEH GYSPFWQSLG KRFFSMDFSR ADFLCGTGQK AFIAELMPKH PIYTHFLSQE AQDVIQVHP QTAPARAVLE KEGFRYRNYI DIFDGGPTLE CDIDRVRAIR KSRLVEVAEG QPAQGDFFAC LVANENYHHF RVVLRVTDPA TERLILTAAQ LDALKCHAGD RVRLVRLCAE EKTA
<b>Source</b>	Baculovirus
<b>Target Names</b>	astA
<b>Protein Names</b>	Recommended name: Arginine N-succinyltransferase Short name= AST EC= 2.3.1.109 Alternative name(s): AOST
<b>Expression Region</b>	1-344
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.