



Recombinant Escherichia coli O17:K52:H18 L-arabinose isomerase (araA)

Product Code	CSB-MP484574EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7N7T6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTIFDNYEVW FVIGSQHLYG PETLRQVTQH AEHVVNALNT EAKLPCKLVL KPLGTTTPEI TAICRDANYD DRCAGLVVWL HTFSPAKMWI NGLTMLNKPL LQFHTQFNAA LPWDSIDMDF MNLNQTAHGG REFGFIGARM RQQHAVVTGH WQDKQAHERI GSWMRQAVSK QDTRHLKVCR FGDNMREVAV TDGDKVAAQI KFGFSINTWA VGDLVQVVNS ISDGDVNALV DEYESCYTMT PATQIHGEKR QNVLEAARIE LGMKRFLEQG GFHAFTTTTFE DLHGLKQLPG LAVQRLMQQG YGFAGEGDWK TAALLRIMKV MSTGLQGGTS FMEDYTYHFE KGNDLVLGSH MLEVCPSSIAA EEKPILDVQH LGIGGKDDPA RLIFNTQTGP AIVASLIDLG DRYRLLVNCI DTVKTPHSLP KLPVANALWK AQPDLPTASE AWILAGGAHH TVFSHALNLN DMRQFAEMHD IEITVIDNDT RLPFAKDALAR WNEVYYGFRR
Source	Mammalian cell
Target Names	araA
Protein Names	Recommended name: L-arabinose isomerase EC= 5.3.1.4
Expression Region	1-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.