



# Recombinant Escherichia coli O17:K52:H18 Cysteine--tRNA ligase (cysS)

<b>Product Code</b>	CSB-MP484586EOG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7N980
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLKIFNTLTR QKEEFKPIHA GEVGMVCGI TVYDLCHIGH GRTFVAFDVV ARYLRFLGYK LKYVRNITDI DDKIHKRANE NGESFVALVD RMIAEMHKDF DALNILRPDM EPRATHHIAE IIELTEQLIA KGHAYVADNG DVMFDVPTDP TYGVLSRQDL DQLQAGARVD VVDDKRNPMDFVLWKMSKEG EPSWSPWGP GRPGWHIECS AMNCKQLGNH FDIHGGGSDL MFPHHENEIA QSTCAHDGQY VNYWMHSGMV MVDREKMSKS LGNFFTVRDV LKYYDAETVR YFLMSGHYRS QLNyseenlk QARAALERLY TALRGTDKTV APAGGEAFEA RFIEAMDDDF NTPEAYSVLF DMAREVNRLK AEDMAAANAM ASHLRKLSSV LGLLEQEPEA FLQSGAQADD SEVAEIEALI QQRLDARKAK DWAAADAARD RLNEMGIVLE DGPQGTTWRR K
<b>Source</b>	Mammalian cell
<b>Target Names</b>	cysS
<b>Protein Names</b>	Recommended name: Cysteine--tRNA ligase EC= 6.1.1.16 Alternative name(s): Cysteinyl-tRNA synthetase Short name= CysRS
<b>Expression Region</b>	1-461
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.