



# Recombinant Escherichia coli O7:K1 Ribosomal RNA large subunit methyltransferase F (rlmF)

<b>Product Code</b>	CSB-MP484686EON
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7NNN9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSAQKPG LHP RNRHHSRYDL ATLCQVNP EL RQFLTLTPAG EQSVDFANPL AVKALNKALL THFYAVANWD IPDGFLCPPV PGRADYIHYL ADLLAEATGS IPANASILDI GVGANCIYPL IGVHEYGWRF TGSETSSQAL SSAQAIISAN PGLNRAIRLR RQKESGAIFN GIIHKNEQYD ATLCNPPFHD SAAAARAGSE RKRRNLGLNK DDALNFGGQQ QELWCEGGEV AFIKKMIEES KGFAKQVMWF TSLVSRGENL PPLYRALTDV GAVKVVKKEM AQQQKQSRFI AWTFMNDQQR RRFVNRQR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rlmF
<b>Protein Names</b>	Recommended name: Ribosomal RNA large subunit methyltransferase F EC= 2.1.1.181 Alternative name(s): 23S rRNA mA1618 methyltransferase rRNA adenine N-6-methyltransferase
<b>Expression Region</b>	1-308
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.