



# Recombinant Escherichia coli O7:K1 NADH-quinone oxidoreductase subunit C/D (nuoC)

<b>Product Code</b>	CSB-EP484689EON-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7NNW5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVNNMTDLTA QEPAWQTRDH LDDPVIGELR NRFGPDAFTV QATRTGVPVV WIKREQLLEV GDFLKKLPKP YVMLFDLHGM DERLRTHREG LPAADFSVFY HLISDRNRD IMLKVALAEN DLHVPTFTKL FPNANWYERE TWDLFGITFD GHPNLRRIMM PQTWKGHPLR KDYPARATEF SPFELTKAKQ DLEMEALTFK PEEWGMKRGT ENEDFMFLNL GPNHPSAHGA FRIVLQLDGE EIVDCVPDIG YHHRGAEKMG ERQSWHSYIP YTDRIEYLG GCVNEMPYVLA VEKLAGITVP DRVNVIRVML SELFRINSHL LYISTFIQDV GAMTPVFFAF TDRQKIYDLV EAITGFRMHP AWFRI GGVAH DLPRGWDRLL REFLDWMPKR LASYEKAALQ NTILKGRSQG VAA YGAKEAL EWGTTGAGLR ATGIDFDVRK ARPYSGYENF DFEIPVGGGV SDCYTRVMLK VEELRQSLRI LEQCLNNMPE GPFKADHPLT TPPPKERTLQ HIETLITHFL QVSWGPVMPA NESFQMIEAT KGINSYYLTS DGSTMSYRTR IRTPSY AHLQ QIPAAIRGSL VSDLIVYLG S IDFVMSDVDR
<b>Source</b>	E.coli
<b>Target Names</b>	nuoC
<b>Protein Names</b>	Recommended name: NADH-quinone oxidoreductase subunit C/D EC= 1.6.99.5 Alternative name(s): NADH dehydrogenase I subunit C/D NDH-1 subunit C/D
<b>Expression Region</b>	1-600
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.