



Recombinant Escherichia coli O127:H6 Probable crotonobetaine/carnitine-CoA ligase (caiC)

Product Code	CSB-BP484757EOB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7UI83
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Purity	>85% (SDS-PAGE)
Sequence	MDIIGGQHLR QMWDDLADV GHKTALICES SSGVVNRYSY LELNQEINRT ANLFYTLGIR KGDKVALHLD NCPEFIFCWF GLAKIGAIMV PINARLLREE STWILQNSQA CLLVTSAQFY PMYQQIQQED ATQLRHICLT DVSLPADDGV SSFTQLKNQQ PATLCYAPPL STDDTAEILF TSGTTSRPKG VVITHYNLRF AGYSAWQCA LRDDDVYLTV MPAFHIDCQC TAAMAAFSTG ATFVLVEKYS ARAFWGQVQK YRATITECIP MMIRTLMVQP LSANDQQHRL REVMFYLNLS EQEKDAFCER FGVRLLSYG MTETIVGIIG DRPGDKRRWP SIGQAGFCYE AEIRDDHNRP LPAGEIGEIC IKGVPGKTIF KEYFLNPKAT AKVLEADGWL HTGDTGYRDE EGFFYFVDRR CNMIKRGGEN VSCVELENII ATHPKIQDIV VVGIKDSIRD EAIKAFVVLN EGETLSEEEF FCFCEQNMAK FKVLSYLEIR KDLPRNCSGK IIRKNLK
Source	Baculovirus
Target Names	caiC
Protein Names	Recommended name: Probable crotonobetaine/carnitine-CoA ligase EC= 6.2.1.-
Expression Region	1-517
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.